VETLIS® Brucella iELISA Ovinos

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VETLIS®
Brucella iELISA Ovines

Presentations / Components

| Components | Description | 15-VL04-2P 2 microplates | 15-VL04-5P 5 microplates |
|------------------------------|---|---|---------------------------------------|
| Microplate | 96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigen, sealed and stored in dry | 2 microplates (192 tests) | 5 microplates (480 tests) |
| Conjugate | Lyophilized (horseradish peroxidase conjugated anti-bovine IgG antibodies) | 2 x 13 ml 4 x 10 ml | 5 x 13 ml 4 x 10 ml |
| Substrate-chromogen solution | Ready to use [3,3',5,5'-Tetramethyl-bencidine (H ₂ O ₂)] – STORE IN THE DARK | 1 x 30 ml | 1 x 60 ml |
| Stop solution | Ready to use (contains HCl 1%) - CORROSIVE | 1 x 30 ml | 1 x 60 ml |
| Wash solution | 10X concentrate | 1 x 120 ml, 10X (c. s. p. 12 litros) | 2 x 120, 10X (c. s. p. 2,4 litros) |
| Positive control | Serum - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%] | 1 x 0.1 ml | 2 x 0.2 ml |
| Negative control | Serum - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%] | 1 x 0.1 ml | 2 x 0.2 ml |
| Polyethylene bag | With ziplock and reusable | 1 | 1 |
| Instruction manual | | 1 | 1 |

Name and application

VETLIS®
Brucella iELISA Ovine

Indirect ELISA kit for the detection of anti *Brucella ovis* antibodies in ovine serum.

The kit VETLIS® Brucella iELISA Ovine is an ELISA (Enzyme Linked Immunosorbent Assay) kit for the indirect detection of specific antibodies against the rough lipopolysaccharide (rLPS) of *Brucella ovis* in ovine serum. Thanks to its exclusive Pure-R technology, the kit has an excellent diagnostic performance. VETLIS® Brucella iELISA Ovinos can be used as a screening and/or a confirmatory test.

Introduction

Brucella ovis is the causative agent of "Ovine Epididymitis", a transmissible infectious disease that affects only sheep. In ovine the disease occurs clinically or subclinically and is characterized by genital lesions in males and placentitis in females. Therefore, the main consequences of the disease are a decrease fertility in rams, abortions in sheep (rare) and increase perinatal mortality. In particular, the control of ovine epididymitis depends exclusively on the detection and killing of infected animals. The confirmatory diagnosis of brucellosis is made directly by isolating the microorganism from samples of semen or tissues of the rams, or vaginal fluids and milk from the sheep. However, the slow growth of *Brucella* from primary cultures (up to 7 days), the risk involved in its management and the low sensitivity of the bacteriological isolation determine that the diagnosis based exclusively on the isolation of *Brucella* is not always feasible and effective. For these reasons, the laboratory diagnosis of the disease is based mainly on serological diagnosis by detecting specific antibodies against the infectious agent in serum samples.

Principles of the technique

VETLIS® Brucella ELISA Ovine kit is a solid phase enzyme immunoassay for the detection of IgG specific antibodies against the rough lipopolysaccharide (rLPS) of *Brucella ovis* in serum samples. In this assay the antibodies present in the sample react with the LPSr (purified by the exclusive technology Pure-R) that coat the wells. If the sample contains antibodies against *B. ovis* LPSr, they bind to the well. Upon addition of the horseradish peroxidase (HRP) conjugated antibody, a complex is formed with the IgG antibodies attached to the antigen. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate and color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigens. The higher the antibody titer present in the sample, the more intense color develops in the test wells. Adding the stop solution, upon which a color change to yellow is observed, stops the reaction. The results are measured in a microplate spectrophotometer determining the absorbance (Abs) at 450 nm.

Storage and expiration

Store: between 2 and 8°C.

Transport temperature: 4 to 15°C for 72 hours. Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

- High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).
- Containers for multichannel pipette.
- Disposable pipette tips.
- Microplate reader for absorbance determination at 450 nm.

- Orbital rotator.
- Timer.
- Tubes for the dilutions of the controls and samples.
- Container/graduated cylinder for the preparation of the wash solution.
- Distilled (Type II) or ultrapure (Type I) water.
- Microplate covers (lid, aluminium foil or adhesive).

Precautions

1. Carefully read and strictly follow all instructions detailed in the Instruction manual included in the kit.
2. Store the kit and all reagents at 2-8 °C.
3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.
4. Handle all materials and reagents according to the Good Laboratory Practices.
5. Serum and milk samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.
6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.
7. Do not mix components or instruction manuals from different test kits batches.
8. Care should be taken to avoid contamination of kit components.
9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.
10. Use different pipettes tips for each re-

agent and for each sample.

11. Include positive and negative controls, and the target sample in duplicate in each test.
12. For reconstitution of the reagents use only distilled (Type II) or ultrapure (Type I) water.

Reagents preparation

a. Wash solution 1X

The 10X concentrated Wash solution should be brought to room temperature (20-25°C) and mixed to ensure dissolution of any precipitated salts.

The 10X concentrated Wash solution must be diluted 1/10 with distilled (Type II) or ultrapure (Type I) water; for example, adding 30 ml of 10X Wash solution to 270 ml of water (sufficient quantity of Wash solution 1X to process one microplate: 13 ml to reconstitute the IgG conjugate and the rest for the washes and preparation of the controls and samples). Mix very well before using. Once prepared, the 1X Wash solution can be stored at 2-8 °C for up to 7 days.

It is recommended to prepare sufficient quantity to process the plates of the day.

b. Conjugate

Option 1: when it is required to process the entire plate, reconstitute the corresponding lyophilized conjugate (13 ml vial) by adding 13 ml of the previously prepared 1X Wash Solution (see previous item).

For greater accuracy, it is recommended to pipet the 13 ml using two graduation marks of the pipette (for example: pipetting twice from zero to 6.5 ml).

Option 2: when it is not required to process the entire plate, reconstitute the corresponding lyophilized conjugate (10 ml vial) by adding 10 ml of the previously prepared 1X Wash Solution (see previous item).

For greater accuracy, it is recommended to pipet the 10 ml using two graduation marks of the pipette (for example: pipetting twice from zero to 5.0 ml).

-Reconstitute the conjugate only immediately before use.

-Add the 1X Wash Solution carefully. Let the bottle rest for one minute and mix very well (cover the bottle and mix twice by shaking head-tail).

-The remaining conjugate, once reconstituted, must be discarded.

Preparation of controls

a. Positive and negative controls

Controls should be diluted 1:200 using the Wash solution 1X.

In an Eppendor-type tube add 5 μ l of the corresponding control to 1000 μ l of Wash solution 1X (dilution 1:200) and mix well. Use 100 μ l of each diluted control per well.

Note: controls should be analyzed at least in duplicate.

b. Sample blank (SB)

On each plate a SB consisting of Wash solution 1X should be added. Dispense directly and in duplicate 100 μ l per well of Wash solution 1X.

Sample preparation

Serum samples should be diluted 1:200 using the Wash solution 1X.

In an Eppendor-type tube add 5 μ l of serum to 1000 μ l of Wash solution 1X (dilution 1:200) and mix well. Use 100 μ l of the diluted sample per well.

Notes:

-Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in a refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bringing them to room temperature, and homogenized before performing the analysis.

-It is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

Test procedure

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature. If the entire plate is not used, separate only the 8-well strips needed to analyze the controls and samples. Store the rest of the strips, together with the desiccants, in the polyethylene bag with hermetic seal included in the kit and store it again in the refrigerator (2-8°C).

Note: once the original packaging of the plate has been opened, the strips can be kept in the closed polyethylene bag and in the refrigerator (2-8 °C) up to 90 days.

2. Add 100 μ l per well and in duplicate the positive and negative controls and the sample blank. See Item Preparation of controls.

3. Add 100 μ l per well of each sample. See Item Sample preparation.

Note: samples can be processed in a single well or in duplicate. However, to obtain more reliable results it is recommended to run samples at least in duplicates.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 4 times with 200 μ l of 1X Wash Solution per well. See item Reagents preparation.

Note: remove the liquid completely after each washing by dump and tapping onto absorbent material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 μ l of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 μ l of Substrate-chromogen solution to each well, and incubated with gentle orbital shaking for 10 minutes (\pm 1

min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

Important: respect incubation times.

Calculations

Expression of results

Results are expressed as the percentage of reactivity ($PR_{PC} = 100\%$) with respect to the positive control (PC) included in each test. The absorbance values measured at 450 nm (Abs_{450}) of each sample are related to the Abs_{450} value of the PC as follows:

$$PR_{SAMPLE} = \frac{Abs_{450} \text{ Sample}^*}{Average Abs_{450} \text{ PC}} \times 100$$

* If the sample was run in duplicate, consider the average of the Abs_{450} value.

Validity criteria

To confirm the validity of the test the following criteria must be met:

- The average value of the PC (Abs_{450}) must be higher than 1.1 ($Abs_{450} \text{ PC} > 1.1$).
- The Abs_{450} values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.
- The Abs_{450} values of the sample blank must be less than 0.11 ($Abs_{450} \text{ SB} < 0.11$).

- The average Abs_{450} value of the NC must be less than 0,15 (Average $Abs_{450} \text{ NC} < 0.15$).

- The Average $Abs_{450} \text{ PC} / \text{Average } Abs_{450} \text{ NC}$ ratio must be higher than 7 (Average $Abs_{450} \text{ PC} / \text{Average } Abs_{450} \text{ NC} > 7$).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the Instruction manual included in the kit.

Interpretation of results

The test sample results should be interpreted as follows:

Serum (dilution 1:200)

| PR | Interpretation |
|--------------|----------------|
| ≤ 19,5% | No Reactive |
| > 19,5 ≤ 50% | Indeterminate |
| > 50% | Reactive |

In the case of obtaining an "Indeterminate" result, the sample should be re-tested. If the result is still indeterminate, it is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is indeterminate again, the epidemiological situation should be considered and the sample should be analyzed with another technique if available.

To facilitate the calculations, the evaluation of the validity criteria of the test and the interpretation of the results, a calculation sheet is available on the company's website (<https://www.chemtest.net/elisa-ovinos.php>) or request it via email to info@chemtest.net. **Important:** always check that the version number of the instruction manual provided in the kit matches the version of the spreadsheet.

Sensitivity and specificity

More than 350 serum samples obtained from animals infected with *B. ovis* and

from healthy animals were analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay. Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed. This analysis allowed us to select the cut-off value that concurrently optimized the Se and Sp (19.5%) of the assay and the cut-off values for which a 100 % Se or Sp was achieved (see Table 1). The cut-off value that simultaneously maximized the Se and Sp (maximum value of the Youden's index, J) matches the cut-off value for which the Se is 100%.

Table 1. Sensitivity, specificity and Youden's index.

| Cut-off (PR) ^a | Se (%) ^b | Sp (%) ^b | J ^c |
|---------------------------|-----------------------|-----------------------|----------------|
| > 19,5 | 100 (93,5 - 100) | 90,2 (84,3 - 94,4) | 0.902 |
| > 50 | 83,6 (71,2 - 92,3) | 100 (97,6 - 100) | 0.836 |

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Chemtest Argentina S.A.