# VETLIS® Brucella Glyco-iELISA Bovinos

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**VETLIS®**  
**Brucella Glyco-iELISA Bovines (serum/milk)**

**Presentations / Components**

| Components                   | Description                                                                                                        | 15-VL01-2P<br>2 microplates              | 15-VL01-5P<br>5 microplates         |
|------------------------------|--------------------------------------------------------------------------------------------------------------------|------------------------------------------|-------------------------------------|
| Microplate                   | 96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigen, sealed and stored in dry  | 2 microplates<br>(192 tests)             | 5 microplates<br>(480 tests)        |
| Conjugate 20X IgG            | 20X concentrated (horseradish peroxidase conjugated anti-bovine IgG antibodies)                                    | 1 x 1.25 ml                              | 1 x 3.2 ml                          |
| Conjugate diluent            | Ready to use - CONTAINS PRE-SERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> Hg-NaO <sub>2</sub> S) 0.01%].    | 1 x 25 ml                                | 1 x 60 ml                           |
| Substrate-chromogen solution | Ready to use [3,3',5,5'-Tetramethyl-benzidine (H <sub>2</sub> O <sub>2</sub> )] – STORE IN THE DARK                | 1 x 25 ml                                | 1 x 60 ml                           |
| Stop solution                | Ready to use (contains HCl 1%) - CORROSIVE                                                                         | 1 x 25 ml                                | 1 x 60 ml                           |
| Wash solution                | 10X concentrate                                                                                                    | 1 x 120 ml, 10X<br>(c. s. p. 1.2 liters) | 1 x 120 ml<br>(c. s. p. 1.2 liters) |
| Sample diluent               | To reconstitute - CONTAINS PRE-SERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> Hg-NaO <sub>2</sub> S) 0.01%]. | 1 x 120 ml                               | 2 x 120 ml                          |
| Positive control             | Serum - CONTAINS PRESERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> HgNaO <sub>2</sub> S) 0.01%]              | 1 x 0.1 ml                               | 1 x 0.2 ml                          |
| Negative control             | Serum - CONTAINS PRESERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> HgNaO <sub>2</sub> S) 0.01%]              | 1 x 0.1 ml                               | 1 x 0.2 ml                          |
| Polyethylene bag             | With ziploc. Reusable                                                                                              | 1                                        | 1                                   |
| Instruction manual           |                                                                                                                    | 1                                        | 1                                   |

## Name and application

### VETLIS® Brucella Glyco-iELISA Bovines (serum/milk)

Indirect ELISA kit for the detection of anti Brucella antibodies in bovine serum and milk.

The kit VETLIS® Brucella Glyco-iELISA Bovines is an iELISA kit (Enzyme Linked Immunosorbent Assay) is exclusively based in the detection of specific antibodies against the O polysaccharide of the lipopolysaccharide (LPS) of *Brucella abortus*, *Brucella melitensis* and *Brucella suis* in bovine serum and milk.

Thanks to its exclusive GlycoEng technology, VETLIS® Brucella Glyco-iELISA Bovines has an excellent diagnostic performance and can be used as a screening or confirmatory test.

Anti-O polysaccharide antibodies are an excellent specific marker of brucella infection. The exclusive detection of these antibodies allows to confirm the diagnosis of infection by eliminating the risk of false-positive cross-reactions due to the presence of antibodies directed against common epitopes present in the lipid A and the core of the LPS from other bacteria. In addition, false positive reactions due to the presence of residual antibodies produced in response to vaccination with S19 are minimized. In some cases, samples obtained before 6 months post-vaccination with the strain S19 strain may be reactive. Therefore, the samples from animals vaccinated with S19 should be obtained and analyzed after 6 months post-vaccination.

VETLIS® Brucella Glyco-iELISA Bovines is a DIVA (Differentiating Infected from Vaccinated Animals) test for rough vaccine strains such as RB51 and Δpgm because it exclusively detects anti-O polysaccharide antibodies, an antigen that is absent in these vaccine strains.

## Introduction

Brucellosis is a highly contagious zoonotic disease caused by Gram-negative bacteria of the genus *Brucella*. The species *Brucella abortus* is the main causative agent of brucellosis in cattle. This disease causes abortions in advanced gestation periods in pregnant females, infertility, placental retention, premature births and decreased milk production. In males, the infection can cause infertility, epididymitis and orchitis. Additionally, detection of the disease in a region or country imposes restrictions on the movement and sale of livestock. For these reasons, bovine brucellosis causes great economic losses and, therefore, should be monitored in herds to detect infection. Due to the high economic impact of the disease in animal production and health, and to the risk of transmission to humans, most countries have implemented programs to control and/or eradicate brucellosis in cattle. Brucellosis in humans can be debilitating and disabling and is an important public health problem. In the absence of a vaccine against human brucellosis, disease prevention depends primarily on the control of brucellosis in animals, the natural reservoir of the disease. In particular, the control of bovine brucellosis depends primarily on the vaccination with vaccine strain S19, in combination with the detection and slaughter of infected animals. The confirmatory diagnosis of brucellosis is performed directly by isolating the organism from blood cultures, bone marrow or other tissues. However, due to the slow growth of *Brucella* in primary cultures (up to 7 days), the risk involved in manipulating the bacteria and low sensitivity of its isolation, determined that the isolation of *Brucella* for diagnosis is not always feasible and effective. For this reason, the laboratory diagnosis of bovine brucellosis is mainly based on the detection of specific antibodies against the infectious agent in serum and/or milk samples. In this sense, it is essential

that the serological method used allows discrimination of vaccinated from naturally infected animals, and thus avoid unnecessary slaughter of healthy animals.

## Principles of the technique

The VETLIS® Brucella Glyco-iELISA Bovines kit is the first indirect solid phase enzyme immunoassay based exclusively on the detection of *Brucella abortus*, *Brucella melitensis* and *Brucella suis* anti-O polysaccharide IgG antibodies. A positive result by Brucella Glyco-iELISA Bovines indicates the presence of anti-O polysaccharide antibodies.

In this assay, serum or milk samples are exposed to the wells coated with the Brucella antigen (purified recombinant glycoprotein). If the sample contains anti *Brucella abortus* O polysaccharide antibodies, they bind to the antigen in the well. Upon addition of the conjugate, horseradish peroxidase (HRP) conjugated anti-bovine IgG antibodies, a complex is formed with the bovine IgG antibodies attached to the antigen in the well. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate in those samples that are positive. The color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigen. The higher the antibody titer present in the sample, the more intense color develops in the test wells. The reaction is stopped by adding the stop solution, upon which a color change to yellow is observed. The reading of the results is performed in a microplate spectrophotometer measuring the absorbance (Abs) at 450 nm.

## Storage and expiration

Store: between 2 and 8°C.

Transport temperature: 4 to 15°C for 72 hours. Proper storage of the kit and all

reagents will allow the use until the expiration date.

Expiration: 12 months. Do not use reagents from different lots or after the expiration date.

## Materials needed but not provided

- High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).
- Containers for multichannel pipette.
- Disposable pipette tips.
- Microplate reader for absorbance determination at 450 nm.
- Orbital rotator.
- Timer.
- Tubes for the dilution of the controls, samples and conjugate.
- Container (graduated cylinder) for the preparation of the wash solution.
- Distilled (Type II) or ultrapure (Type I) water or any similar high quality water.
- Microplate covers (lid, aluminium foil or adhesive).

## Precautions

1. Carefully read and strictly follow all instructions detailed in the instruction manual included in the kit.
2. Store the kit and all reagents at 2-8 °C.
3. All reagents should equilibrate at room temperature (20-25°C) before use and return to 2-8 °C following use.
4. Handle all materials and reagents according to the Good Laboratory Practices.
5. Serum and milk samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.
6. Do not use the kit beyond expiration date or if the kit components were not

kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the blank sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled (Type II) or ultrapure (Type I) water.

## Reagents preparation

### a. Wash solution 1X

The 10X concentrated Wash solution should be brought to room temperature (20-25°C) and mixed to ensure dissolution of any precipitated salts.

The 10X concentrated Wash solution must be diluted 1/10 with distilled (Type II) or ultrapure (Type I) water; for example, adding 30 ml of 10X Wash solution to 270 ml of water (sufficient quantity of Wash solution 1X to process one microplate: 120 ml to reconstitute the Sample diluent and the rest for the washes). Mix very well before using. Once prepared, the 1X Wash solution can be stored at 2-8 °C for up to 7 days.

Note: If crystals appear in the 10X Wash solution, the bottle must be heated not above 45°C and shake well until completely dissolve.

It is recommended to prepare sufficient quantity to process the plates of the day.

### b. Conjugate

Dilute 1/20 using the Conjugate diluent. For example, to process the entire plate add 9.5 ml of the Conjugate diluent and 0.5 ml of the 20X IgG conjugate. Homogenize very well.

#### Important:

-Do not remove the 20X IgG Conjugate from the refrigerator until use.

-Shake the Conjugate diluent thoroughly before use.

-Perform the dilution of the conjugate immediately before step 5 (see item Test procedure). The Conjugate diluent should be at room temperature.

-Only prepare enough volume of diluted conjugate to process the required number of wells.

-The remaining solution of the diluted conjugate should be discarded.

### c. Sample diluent

Reconstitute the Sample diluent by adding the 1X Wash solution previously prepared (see previous item) up to the level indicated in the corresponding bottle (approximately 120 ml) and mix very well. Prepare immediately before use.

The remaining Sample diluent can be stored in the refrigerator (2 to 8°C) up to 3 months or at -20°C for a year. Thaw it completely, bring to room temperature and mix very well until complete dissolution before using.

## Preparation of controls

### a. Positive (PC) and negative (NC) controls

Controls (PC and NC) should be diluted 1:200 using the Sample diluent.

Add 2 µl of the corresponding control to 400 µl of Sample diluent (dilution 1:200) and mix well. Use 100 µl of each diluted control per well and in duplicate.

Notes:

-Shake the Sample diluent thoroughly before preparing the controls.

-The provided positive and negative controls are used to process either serum or milk samples.

### b. Sample blank (SB)

On each plate, a sample blank consisting of Sample diluent should be added. Dispense directly 100 µl per well of Sample diluent.

## Sample preparation

### a. Serum samples

Serum samples should be diluted 1:200 using the Sample diluent. Add 2 µl of serum to 400 µl of Sample diluent. Use 100 µl of diluted sample per well.

Notes:

- Shake the Sample diluent thoroughly before preparing the samples.

- Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bring to room temperature and homogenize before performing the analysis.

### b. Milk samples

Individual or tank milk samples (pooled samples from the storage tank) are processed undiluted. 100 µl of sample is used per well. Freshly collected milk samples, samples stored in a refrigerator (2-8°C) for no more than 3 days, or frozen milk samples (-20°C) can be used. If samples stored in a refrigerator or frozen are used, they must be brought to room temperature and thoroughly homogenized before analysis.

General Notes: Samples (serum, individual milk or tank milk) can be analyzed in a single well. However, it is recommended to analyze the samples at least in duplicate. Samples contaminated and/or not well conserved may produce erroneous results.

## Test procedure

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature. If the entire plate is not used, separate only the 8-well strips needed to analyze the controls and samples. Store the rest of the strips, along with the desiccants, in the original foil bag (seal the open side with a couple of folds and adhesive tape) and place it in the ziplock-type

polythene bag included in the kit and store in the refrigerator (2-8°C).

Note: once the original packaging of the plate has been opened, the strips can be kept in the closed polyethylene bag and in the refrigerator (2-8 °C) up to 90 days.

2. Add in duplicate 100 µl per well of the positive control, negative control and the blank sample. See Item Preparation of controls.

3. Add 100 µl per well of each serum or milk sample. See Item Sample preparation.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 4 times with 200 µl of 1X Wash Solution per well. See item Reagents preparation.

Note: Remove the liquid completely after each washing by dump and tapping onto absorbent material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 µl of Substrate-chromogen solution to each well, and incubated with gentle orbital shaking for 10 minutes ( $\pm 1$  min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can

lead to the appearance of a black precipitate due to the precipitation of the substrate.

Important: respect incubation times.

Calculations

Expression of results

Results are expressed as the percentage of reactivity (PR) with respect to the positive control included in each test. The absorbance values measured at 450 nm (Abs450) of each sample are related to the Abs450 value of the positive control (PC) as follows:

$$PR_{SAMPLE} = \frac{Abs450\ Sample^*}{Average\ Abs450\ PC} \times 100$$

\* If the sample was run in duplicate, consider the average of the Abs450 value.

Validity criteria

To confirm the validity of the test, the following criteria must be met:

- The average value of the PC (Abs450) must be higher than 1.8 (Abs450 PC > 1.8).
- The Abs450 values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.
- The Abs450 values of the blank sample must be less than 0.2 (Abs450 Blank sample < 0.2).
- The average Abs450 value of the NC must be less than 0,25 (Average Abs450 NC < 0.25).
- The Average Abs450 PC / Average Abs450 NC ratio must be higher than 7,2 (Average Abs450 PC / Average Abs450 NC > 7,2).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the Instruction manual included in the kit.

Interpretation of results

Serum samples (dilution 1:200): the results should be interpreted as follows.

| PR         | Interpretation |
|------------|----------------|
| ≤ 14%      | Negative       |
| > 14 ≤ 24% | Indeterminate  |
| > 24%      | Positive       |

**Positive result:** when the PR is higher than 24%.

**Negative result:** when the PR is lower than or equal to 14%.

**Indeterminate result:** when the PR is higher than 14% and lower or equal to 24%.

Tank milk samples (undiluted): the results should be interpreted as follows.

| PR           | Interpretation |
|--------------|----------------|
| ≤ 25%        | Negative       |
| > 25 ≤ 26.5% | Indeterminate  |
| > 26.5%      | Positive       |

**Positive result:** when the PR is higher than 26.5%.

**Negative result:** when the PR is lower than or equal to 25%.

**Indeterminate result:** when the PR is higher than 25% and lower or equal to 26.5%.

Individual milk samples (undiluted): the results should be interpreted as follows.

| PR         | Interpretation |
|------------|----------------|
| ≤ 45%      | Negative       |
| > 45 ≤ 55% | Indeterminate  |
| > 55%      | Positive       |

**Positive result:** when the PR is higher than 55%.

**Negative result:** when the PR is lower than or equal to 45%.

**Indeterminate result:** when the PR is higher than 45% and lower or equal to 55%.

In the case of obtaining an indeterminate result, the sample should be re-tested. If the result is still indeterminate, it

is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is indeterminate again, the epidemiological situation should be considered and the sample should be analyzed with another technique if available.

To facilitate the calculations, the evaluation of the validity criteria of the test and the interpretation of the results, a calculation sheet is available on the company's website (<https://www.chemtest.net/en/elisa-bovines.php>) or request it via email to [info@chemtest.net](mailto:info@chemtest.net).

**Important:** always check that the version number of the instruction manual provided in the kit matches the version of the spreadsheet.

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#### **Referencias / References**

-Ciocchini, AE; Rey Serantes, DA; Melili, LJ; Guidolin, LS; Iwashkiw, J; Elena, S; Franco, C; Nicola, AM; Feldman, MF; Comerçi, DJ and Ugalde, JE. A bacterial engineered glycoprotein as a novel antigen for diagnosis of bovine brucellosis. *Veterinary Microbiology*, Aug. 2014, 27; 172(3-4):455-65.

-Iwashkiw, J., Fentabil, M., Faridmoayer, A., Mills, D.C., Peppler, M. Czibener, C., Ciocchini, A.E., Comerçi, D.J., Ugalde, J.E. and Feldman, M.F. Exploiting the *Campylobacter jejuni* protein glycosylation system for glycoengineering vaccines and diagnostic tools directed against brucellosis. *Microbial Cell Factories* 2012 Jan 25;11-13.

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