CHEMLIS® Chagas R-iELISA

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CHEMLIS®
Chagas R-iELISA

Kit components / Presentations

Components	Description	REF: 15-CL01-2P Presentation 2 microplates	REF: 15-CL01-4P Presentation 4 microplates
Microplate	96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigens, sealed and stored dry	2 microplates (192 tests)	4 microplates (384 tests)
Conjugate	Lyophilized (horseradish peroxidase conjugated anti-human IgG antibodies)	2 units	4 units
Substrate-chromogen solution	Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hydrogen peroxide (H ₂ O ₂)]- STORE IN THE DARK	1 x 30 ml	1 x 60 ml
Stop solution	Ready to use (contains HCl 1%) - CORROSIVE	1 x 30 ml	1 x 60 ml
Wash solution	10X concentrate	1 x 120 ml, 10X (sufficient for 1.2 liters)	1 x 120 ml, 10X (sufficient for 1.2 liters)
Sample diluent	To reconstitute	1 x 60 ml	1 x 120 ml
Buffer PES	Ready to use. Purple red colour - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%].	1 x 60 ml	1 x 120 ml
Positive control	Inactivated serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.2 ml	2 x 0.2 ml
Negative control	Inactivated serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.2 ml	2 x 0.2 ml
Polyethylene bag	With hermetic closure. Reusable.	1 unit	1 unit
Instruction manual		1 unit	1 unit

Name and application

CHEMLIS® Chagas R-iELISA

Indirect ELISA kit for the detection of antibodies against *Trypanosoma cruzi* in human serum and plasma.

CHEMLIS® Chagas R-iELISA is an indirect ELISA kit (Enzyme Linked Immunosorbent Assay) for the detection of IgG specific antibodies against *Trypanosoma cruzi*, causative agent of Chagas disease, in human serum samples or plasma obtained with heparin, sodium citrate or EDTA. Thanks to its exclusive combination of four recombinant antigens, the kit has an excellent diagnostic performance. CHEMLIS® Chagas R-iELISA can be used as a screening and/or a confirmatory test.

Introduction

American Trypanosomiasis or Chagas disease is a major health and economic problem in Latin America caused by the protozoan parasite *Trypanosoma cruzi*. This parasite is found throughout the American continent in a variety of wild and domestic mammalian reservoirs and transmitted by the bite of infected, hematophagous triatomine insects. Apart from vectorial transmission, humans can become infected with *T. cruzi* by ingestion of contaminated food and fluids, from mother to child during pregnancy and through contaminated blood transfusion or organ transplantation. It is estimated that 10-12 million people are currently infected and that up to 120 million individuals living in endemic areas are at risk of infection, and thus at risk of developing cardiac or gut pathology normally associated with chronic Chagas Disease. Increasing travel and immigration has also brought the risk of *T. cruzi* infection into non-endemic countries such as Spain, Australia and the U.S., where some autochthonous cases of transmis-

sion have already been reported. Diagnosis of Chagas disease is challenging because it is often asymptomatic during the acute phase, and evolves into a chronic stage that has different clinical forms. In addition, and due to a major decline in parasitemia during the chronic phase of the disease, *T. cruzi* detection in blood samples by direct examination, blood-culture or xenodiagnosis is difficult as well as labor- and time-consuming. Therefore, detection of anti-*T. cruzi* antibodies in body fluids (typically blood samples) is still the most effective method for demonstrating direct exposure to the parasite. At present, the most widely used serologic methods are indirect hemagglutination assay (IHA), indirect immunofluorescence (IIF), and enzyme-linked immunosorbent assay (ELISA) using total parasite homogenates or semi-purified antigenic fractions. In spite of their simplicity and low cost, these tests show variations in their reproducibility and reliability that can be attributed to the poor standardization of the antigens. The advent of recombinant DNA and peptide synthesis technologies allowed the production and one-step purification of large amounts of highly pure *T. cruzi* immunodominant antigens. An additional advantage of both recombinant antigens and synthetic peptides is that they minimize specificity problems, one of the major drawbacks of immunodiagnosis of Chagas disease.

Principles of the technique

CHEMLIS® Chagas R-iELISA kit is a solid phase enzyme immunoassay for the detection of IgG antibodies against *Trypanosoma cruzi* in human serum or plasma samples. In this assay the antibodies present in the sample react with the recombinant antigens (exclusive combination of four recombinant proteins) that coat the wells. If the sample contains antibodies against these *Trypanosoma cruzi* proteins, they bind to the well. Upon addition of the horse-

radish peroxidase (HRP) conjugated antibody, a complex is formed with the IgG antibodies attached to the antigen. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate and color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigens. The higher the antibody titer present in the sample, the more intense color develops in the test wells. Adding the stop solution, upon which a color change to yellow is observed, stops the reaction. The results are measured in a microplate spectrophotometer determining the absorbance (Abs) at 450 nm.

Storage and expiration

Storage: between 2 and 8°C.

Transport temperature: 4 to 15°C for 72 hours.

Proper storage of the kit and all reagents will allow the use until the expiration date.

Shelf life: 12 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

- High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).
- Containers for multichannel pipette.
- Disposable pipette tips.
- Microplate reader for absorbance determination at 450 nm.
- Timer.
- Container (graduated cylinder) for the preparation of the wash solution.
- Distilled water (Type II) or ultrapure water (Type I).
- Microplate covers (lid, aluminium foil or adhesive).

Precautions

1. Carefully read and strictly follow all instructions detailed in the leaflet includ-

ed in the kit.

2. Store the kit and all reagents at 2-8 °C.

3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.

4. Handle all materials and reagents according to the Good Laboratory Practices.

5. The samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.

6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the target sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled water (Type II) or ultrapure water (Type I).

Preparation of reagents

a. Wash solution 1X

The 10X concentrated Wash solution should be brought to room temperature (20-25°C) and mixed to ensure dissolution of any precipitated salts.

The 10X concentrated Wash solution must be diluted 1/10 with distilled (Type II) or ultrapure (Type I) water; for example, adding 20 ml of 10X Wash solution to 180 ml of water (sufficient quantity of Wash solution 1X to process one microplate: 13 ml to reconstitute the IgG conjugate and the rest for the washes). Mix very well before using. Once prepared, the 1X Wash solution can be stored at 2-8 °C for up to 7 days.

It is recommended to prepare sufficient quantity to process the plates of the day.

b. Conjugate

Reconstitute the conjugate by adding 13 ml (per flask) of 1X Wash solution previously prepared (see previous item). Add the Washing solution carefully. Allow to stand for a minute and mix well. Prepare immediately before use.

For greater accuracy it is recommended to pipet the 13 ml using two graduation marks of the pipette (for example: pipetting twice from zero to 6.5 ml).

The remaining conjugate can be aliquoted in eppendorf type tubes and stored at -20°C (cannot be stored in the refrigerator) for up to 3 months. During this period of time each aliquot can be thawed and re-frozen one time.

c. Sample diluent

Reconstitute the Sample diluent by adding the total volume of the Buffer PES to the Sample diluent bottle and mix very well. After reconstitution, the Sample Diluent should be light green in color. Prepare immediately before use.

The remaining Sample diluent can be stored in the refrigerator (2 to 8°C) up to 3 months or at -20°C for a year. Thaw it completely, bring to room temperature and mix very well until complete dissolution before using.

Controls

a. Positive (PC) and negative (NC) controls

Add 10 µl of the PC or NC directly to the corresponding well containing 200 µl of Sample diluent (dilution 1:20) and mix completely by pipetting up and down with the micropipette until the color of the Sample Diluent changes from light green to dark green. It is recommended to process the controls in duplicate.

b. Sample blank (SB)

On each plate a sample blank consisting of Sample diluent should be added. Dispense directly and in duplicate 200 µl per well of Sample diluent.

Samples

Add 10 µl of the sample (serum or plasma) directly to the corresponding well containing 200 µl of Sample diluent (dilution 1:20) and mix completely by pipetting up and down with the micropipette until the color of the Sample Diluent changes from light green to dark green.

Notes:

-It is recommended to analyze serum or plasma (obtained using heparin, sodium citrate or EDTA as anticoagulant) samples recently obtained (fresh) and free of turbidity.

-If the samples contain particles or a high content of triglycerides, they should be clarified by centrifugation.

-If the samples are not analyzed immediately, they can be kept refrigerated (2 to 8 °C) for up to 72 hs. For longer storage, store at -20°C (-15 to -20°C), taking into account that during storage and as a consequence of the freezing and thawing process, the antibody titer may change. If samples are frozen, they must be completely thawed, homogenized and centrifuged before testing.

-It is recommended to analyze the samples at least in duplicate.

-Contaminated and/or not well conserved samples may produce erroneous results.

Test procedure

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature. If the entire plate is not used, separate only the 8-well strips needed to analyze the controls and samples. Store the rest of the strips, together with

the desiccants, in the polyethylene bag with hermetic seal included in the kit and store it again in the refrigerator (2-8°C).

Note: once the original packaging of the plate has been opened, the strips can be stored in the hermetically sealed polyethylene bag with the desiccants, and in the refrigerator (2-8°C) for up to 4 months, as long as the indicated expiration date is not exceeded.

2. Add 200 µl per well of the reconstituted Sample diluent. See Item-Reagents preparation.

3. Add in duplicate 10 µl per well of the positive control (PC) and negative control (NC), and mix completely by pipetting up and down with the micropipette until the color of the Sample Diluent changes from light green to dark green. Leave two wells free as sample blank (SB). See Item-Controls.

4. Add 10 µl per well of the samples and mix completely by pipetting up and down with the micropipette until the color of the Sample Diluent changes from light green to dark green. See Item-Samples.

Sample diluent color changes:

	Without sample	Positive control	Negative control	Sample blank	Serum or plasma
Colour	Light green	Dark green	Dark green	Light green	Dark green

5. Cover the plate and incubate at 37 ± 1°C for 30 ± 2 minutes.

6. Wash/rinse the plates 4 times with 200 µl of 1X Wash Solution per well. See item-Reagents preparation.

Note: remove the liquid completely after each washing by dump and tapping onto absorbent material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

7. Add 100 µl of Conjugate to each well. See item-Preparation of the conjugate. Cover the plate and incubate at 37 ± 1°C for 30 ± 2 minutes.

8. Repeat step 6.

9. Add 100 µl of Substrate-chromogen solution to each well, mix thoroughly by lightly tapping the edges of the plate and incubate for 10 ± 1 minute at room tempera-

ture (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

10. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 9).

11. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

Calculations

Expression of results

Results are expressed as the percentage of reactivity ($PR_{PC} = 100\%$) with respect to the positive control included in each test. The absorbance values measured at 450 nm (Abs_{450}) of each sample are related to the Abs_{450} value of the positive control (PC) as follows:

$$PR_{\text{Sample}} = \frac{Abs_{450} \text{ Sample}^*}{\text{Average } Abs_{450} \text{ PC}} \times 100$$

* If the sample was run in duplicate, consider the average of the Abs_{450} value.

Validity criteria

To confirm the validity of the test the following criteria must be met:

- The average value of the PC (Abs_{450}) must be higher than 2.0 ($Abs_{450} \text{ PC} > 2.0$).
- The Abs_{450} values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.
- The Abs_{450} values of the blank sample must be less than 0.25 ($Abs_{450} \text{ Blank sample} < 0.25$).
- The average Abs_{450} value of the NC must

be less than 0.45 (Average $Abs_{450} \text{ NC} < 0.45$).

- The Average $Abs_{450} \text{ PC} / \text{Average } Abs_{450} \text{ NC}$ ratio must be higher than 4.4 (Average $Abs_{450} \text{ PC} / \text{Average } Abs_{450} \text{ NC} > 4.4$).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the Instruction manual included in the kit.

Interpretation of results

The test sample results should be interpreted as follows:

Sample (dilution 1:20)

PR	Interpretation
≤ 50 %	Non-Reactive
> 50 ≤ 60 %	Indeterminate
> 60 %	Reactive

In the case of obtaining an "Indeterminate" result, the sample should be re-tested. If the result is still indeterminate, it is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is indeterminate again, the epidemiological situation should be considered and the sample should be analyzed with another technique if available.

To facilitate the calculations, the evaluation of the validity criteria of the test and the interpretation of the results, a calculation sheet is available on the company's website (<https://www.chemtest.net/chemlis-chagas.php>) or request it via email to info@chemtest.net. **Important:** always check that the version number of the instruction manual provided in the kit matches the version of the spreadsheet.

Performance characteristics

Diagnostic sensitivity and specificity

More than 390 serum samples obtained

from patients with acute or chronic Chagas disease and from healthy donors were analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay. Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed which allowed us to determine the cut-off value for which a diagnostic Se and Sp of 100% were obtained.

Hook effect

Based on the analysis of samples with high titers of specific anti-T. cruzi antibodies, it was determined that the CHEMLIS® Chagas R-iELISA kit does not have a hook effect in the detection of IgG antibodies.

Cross reactivity

Cross-reactivity of CHEMLIS® Chagas R-iELISA was evaluated by testing clinical samples obtained from patients with other conditions. Samples with positive serology obtained from patients with the following diseases or bacterial, parasitic, and viral infections were analyzed:

-Bacterial: Escherichia coli O157, O145 and O121 infections, Human Brucellosis, Tuberculosis, Syphilis, Streptococcus pneumoniae, Klebsiella pneumoniae and Acinetobacter baumannii infections.

-Viral: Dengue, VHB, VHC, HTLV I/II, HIV, Hantavirus.

-Parasitic: Leishmaniasis y Toxoplasmosis.

-Autoimmune diseases: rheumatoid arthritis, Systemic Lupus Erythematosus and autoimmune vasculitis.

CHEMLIS® Chagas R-iELISA has a high level of specificity since no cross-reactivity was observed with any of the samples obtained from patients with other conditions or bacterial, viral or parasitic infections.

Interfering substances

The following potentially interfering substances have no impact on the performance of CHEMLIS® Chagas R-iELISA (final test concentrations of interfering substances are indicated in parenthe-

ses):

- Oxalic acid (0.6 mg/ml)
- Bovine albumin (20 mg/ml)
- Human Immunoglobulin G (hIgG) (2 mg/ml)
- Human Immunoglobulin M (hIgG) (2 mg/ml)
- Human Immunoglobulin A (hIgG) (2 mg/ml)
- Cholesterol (5 mg/ml)
- Hemoglobin (10 mg/ml)
- Bilirubin (2 mg/ml)

Repeatability and Reproducibility

Repeatability and reproducibility of CHEMLIS® Chagas R-iELISA was established using in-house reference panels containing negative specimens and a range of positive specimens. There were no differences observed within-run, between-run, between-lots, between-sites, and between-days.

For technical assistance

Tel: (54-11) 5353-6066
info@chemtest.net
chemtest.net

Referencias / References

-Cortina, María E.; Melli, Luciano J.; Roberti, Mariano; Mass, Mijal I.; Longinotti, Gloria; Tropea, Salvador; Lloret, Paulina; Rey Serantes, Diego; Salomón, Francisco; Lloret, Matías; Caillava, Ana J.; Restuccia, Sabrina; Altcheh, Jaime M.; Buscaglia, Carlos A.; Malatto, Laura; Ugalde, Juan E.; Fraigi, Liliana; Moína, Carlos; Ybarra, Gabriel; Ciocchini, Andrés E. and Comerci, Diego J.
Electrochemical magnetic microbeads-based biosensor for point-of-care serodiagnosis of infectious diseases. *Biosensors & Bioelectronics*; 80; 6-2016; 24-33. <https://doi.org/10.1016/j.bios.2016.01.021>.

Indice de símbolos / Symbol index

	Fecha de vencimiento / Expiration date
	Número de catálogo / Catalog number
	Lote / Batch code
	Para diagnóstico in vitro / For in vitro diagnostic use only
	Límites de temperatura / Temperature limits
	Consultar instrucciones para su uso / See instructions for use
	Número de determinaciones / Number of determinations
	Fabricante / Manufacturer
	Riesgo biológico / Biological risk

 **Chemtest Argentina S. A.**
Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.
info@chemtest.net - chemtest.net

ELISA

Humanos

 Chemtest Argentina S. A.
Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.
info@chemtest.net - chemtest.net