CHEMTEST™

CHEMSTRIP® E. coli O157 / O145

Card Test

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Card Test

Name and Intended use

CHEMSTRIP® *E. coli* O157/O145 is a one-step immunochromatographic test for the rapid and early detection of IgM antibodies against the O polysaccharide of the lipopolysaccharide (LPS) of *Escherichia coli* O157 and O145 in human serum, plasma, and whole blood samples.

Thanks to the incorporation of the exclusive GlycoEng technology that allows to obtain recombinant glycoprotein-based antigens, CHEMSTRIP® *E. coli* O157/O145 has an excellent diagnostic performance, minimizes cross-reaction with other Gram-negative bacteria and it is the first serospecific assay for O157 and O145 that allows to differentiate O157 and O145 *E. coli* infections from those caused by other *E. coli* serogroups.

Introduction

Shiga toxin-producing *Escherichia coli* (STEC) is an important foodborne pathogen associated with sporadic cases and outbreaks of diarrhea, bloody diarrhea (BD) and hemolytic uremic syndrome (HUS), a disorder characterized by microangiopathic hemolytic anemia, thrombocytopenia, and acute renal failure. Diarrhea associated HUS is caused mainly by particular STEC serotypes. STEC strains are characterized by the production of Shiga toxins, Shiga 1 (Stx1) and/or Shiga 2 (Stx2), and the production of these toxins is

CHEMTEST™

essential in the pathogenesis of BD and HUS. *E. coli* O157: H7 is the dominant STEC serotype associated with sporadic cases and outbreaks of BD and HUS in different parts of the world.

STEC associated HUS is the most common cause of acute renal failure in children worldwide and accounts for 90% of all cases of HUS in this population. The incidence rate of post-diarrheal HUS varies according to the country and Argentina shows the highest incidence worldwide with 12-14 cases per 100,000 children less than 5 years of age per year. Most cases of HUS in the world are caused by STEC O157, but other serogroups such as O145 can also cause the disease.

The diagnosis of STEC infections is based on the isolation and characterization of STEC strains, and the detection of free fecal Shiga toxin (FFStx). However, cases of HUS with negative STEC isolation and negative FFSTx are frequently observed. It has been demonstrated that STEC is rapidly cleared during the diarrheal prodrome and often ceases before the onset of HUS; for these reasons, it is often not possible to isolate the pathogen or detect its toxins. So, alternative laboratory diagnostic methods have been used such as the detection of antibodies against the LPS of STEC strains. All available serologic tests use inactivated bacteria, bacterial extracts that contain high concentrations of LPS or purified LPS as antigen. Consequently, these assays show high rates of false positive reactions due to the presence of antibodies directed against epitopes present in the lipid A and the core of LPS shared by different STEC strains as well as other bacteria. Bacteriological methods, alone or in combination with tests for the detection FFSTx, provide evidence of STEC infection in only 20 to 50% of children with a clinical diagnosis of HUS. The combination of bacteriological and/or Shiga toxins detection techniques with the serological analysis using the CHEMSTRIP® *E. coli* O157/O145 kit allows to drastically increase the diagnostic efficiency to more than 80% in patients with STEC O157 infections. The exclusive GlycoEng technology incorporated in the CHEMSTRIP® *E. coli* O157/O145 kit confers maximum specificity and sensitivity for the rapid and early detection of anti-O157 and anti-O145 antibodies.

Test principle

The CHEMSTRIP® *E. coli* O157/O145 is an immunochromatographic test based on the detection of IgM antibodies directed against the O polysaccharide of the lipopolysaccharide (LPS) of *E. coli* O157 and O145.

The test is based on the immunological capture of gold nanoparticles coated with goat anti-human IgM antibodies during its passage through a membrane. The IgM isotype immunoglobulins present in the sample react with the gold particles functionalized with the anti-IgM antibodies and the formed complexes migrate by chromatography to the reaction zone. In this area, the antigens for O157 and O145 (recombinant glycoproteins formed by the O157 or O145 polysaccharide linked to a carrier protein) were immobilized in the test lines O157 and O145, respectively, and donkey anti-goat IgG antibodies were immobilized in the control line (CL).

If the sample contains IgM antibodies directed against the *E. coli* O157 and/or *E. coli* O145 polysaccharide, the nanoparticle-antibody complexes are captured by the corresponding antigen immobilized on the membrane which is visualized as a purple red band. Regardless of the presence in the sample of anti-O157 or anti-O145 IgM antibodies, the complexes that were not captured in the test lines continue to migrate and are captured in the CL by the anti-goat IgG antibodies, leading also to the formation of a red purple band. The appearance of the CL indicates that the chromatography has been developed correctly and under conditions that ensure the antigen-antibody reaction.

CHEMTEST™

Kit components

Components	15-CSC17-10C	15-CSC17-20C
Test card (in trilaminate pouch for protection from light and humidity)	x 10	x 20
Sample diluent (plastic dropper)	1 x 2 ml	1 x 4 ml
Graduated and heparinized capillaries	x 10	x 20
Instructions for use	x 1	x 1

Materials needed but not provided

- Timer.
- Micropipette and disposable tips (for serum, plasma, and blood).
- Lancets or needles (for capillary blood).
- Gauze and disinfectant.

Storage and expiration

Store at room temperature between 15 and 25°C.
Transport temperature: 15 to 25°C.
Shelf life: 12 months.

Warnings

- 1- Carefully read and strictly follow all instructions detailed in the manual included in the kit.
- 2- Store the kit and all reagents at room temperature (15-25°C).
- 3- Keep the test card in their original packaging with the silica packets and hermetically sealed to protect them from light and moisture.
- 4- Remove the test card from its pouch only immediately before using it.
- 5- Handle all materials and reagents according to the Good Laboratory Practices.
- 6- Serum, plasma and blood samples should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The rest of the kit components does not represent any potential risk for health and the environment.
- 7- Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.
- 8- Do not mix components or instruction manuals from different test kits batches.
- 9- Proper storage of the kit and all reagents will allow the use until the expiration date.

Collection and storage of samples

Samples that are contaminated and/or not well conserved can produce erroneous results.

Whole blood obtained by lancet or needle puncture

1. Clean the area to be punctured with a gauze or cotton pad with 70% ethanol.
2. Puncture with a sterile lancet or needle. Use a new lancet/needle for each patient.
3. Touch the drop of blood with the tip of the heparinized capillary and collect up the blood to the second black mark (40 µl) indicated on the capillary tube (see figure 1). Avoid the formation of air bubbles.
4. Analyze immediately after collection.

CHEMTEST™


Figure 1: Scheme of the graduated and heparinized capillaries. The red mark indicates that the capillary is coated inside with sodium heparin. Black graduation marks indicate how far the capillary should be loaded to collect 20 µl and 40 µl of blood.

Whole blood obtained by venipuncture

1. Collect the samples by venipuncture in tubes with anticoagulant (heparin or sodium citrate) and mix well to avoid the coagulation of the blood. Do not use EDTA as an anticoagulant.
2. If the samples are not analyzed immediately, they can be stored in a refrigerator (2-8°C) for 1 or 2 days. In these cases, allow the samples to reach room temperature and homogenize well before performing the analysis. Do not freeze whole blood samples.

Plasma

1. Collect blood samples by venipuncture in tubes with anticoagulant (heparin or sodium citrate). Do not use EDTA as an anticoagulant.
2. Obtain the plasma by centrifugation at 3000 rpm for 10 min.
3. If the samples are not analyzed immediately, they can be kept refrigerated (2 to 8 °C) for up to 72 hs. For longer storage, store at -20°C (-15 to -20°C), considering that during storage and because of the freezing and thawing process, the antibody titer may change. If samples are frozen, they must be completely thawed, homogenized, and centrifuged before testing.

Serum

1. Collect the samples by venipuncture in tubes WITHOUT anticoagulant.
2. After blood coagulation (incubate 15 min at 37°C) obtain serum by centrifugation at 3000 rpm for 10 min.
3. If the samples are not analyzed immediately, they can be kept refrigerated (2 to 8 °C) for up to 72 hs. For longer storage, store at -20°C (-15 to -20°C), considering that during storage and because of the freezing and thawing process, the antibody titer may change. If samples are frozen, they must be completely thawed, homogenized, and centrifuged before testing.

Test procedure (see figure 2)

1. Remove the test card from its original packaging and place it on a clean and flat surface.
2. Add the sample and sample diluent to the pad indicated with S.

For whole blood obtained by lancet or needle puncture

-Discharge 10 µl of whole blood onto the test strip by touching the pad (S) with the tip of the preloaded capillary. Wait for the sample to be completely absorbed.

Note: It is recommended to load the capillary to the 40 µl mark and then unload to the middle between the 40 and 20 µl marks (see figure 1).

-Add 3 drops of the Sample diluent on the pad (S) of the reactive strip. Always wait for the volume of diluent added between drops to be absorbed.

For whole blood obtained by venipuncture, serum and plasma

CHEMTEST™

- Add with a micropipette 10 µl of whole blood obtained by venipuncture on the pad (S). Wait for the sample to be completely absorbed.
- Add 3 drops of the Sample Diluent on the pad (S) of the reactive strip. Always wait for absorption of the drop before adding the next one.

3. Wait 10 minutes and read the result.


Figure 2: Instructions for the analysis.

Interpretation and validity criteria of the test (see figure 3)

O157 and O145 positive: when three red purple bands are observed. Two corresponding to the O157 and O145 test lines and the other to the control line (CL). Regardless of the color intensity of the test lines the result should be considered positive.

O157 positive and O145 negative: when two red purple bands are observed. One corresponding to the O157 test line and the other to the control line (CL). Regardless of the color intensity of the test line the result should be considered positive.

O157 negative and O145 positive: when two red purple bands are observed. One corresponding to the O145 test line and the other to the control line (CL). Regardless of the color intensity of the test line the result should be considered positive.

O157 and O145 negative: when only one purple band corresponding to the CL is observed.

Invalid: when no bands are observed or only one band or two bands corresponding to the test lines are observed and the LC is not observed. In this case, the test must be repeated strictly following the procedure detailed in the kit's Instruction Manual.

CHEMTEST™


Figure 3: Interpretation of the results. POS, positive; NEG, negative; INV, invalid.

Important: any line that appears after the indicated reading time, will have no diagnostic value.

Performance evaluation

Evaluation of diagnostic sensitivity (Se) and specificity (Sp)

To determine the diagnostic sensitivity and specificity of the CHEMSTRIP® *E. coli* O157/O145 kit, samples obtained from children under 12 years old were analyzed. The panel of positive samples included samples from patients with clinical diagnosis of acute bloody diarrhea (ABD) or Hemolytic Uremic Syndrome (HUS) and with positive isolation and/or positive serology (assessed using the CHEMLIS® Combi *E. coli* Glyco-iELISA kit) for *E. coli* O157 or O145. The panel of negative samples included samples obtained from healthy children or patients with unrelated pathologies of the same age. Based on these results, the diagnostic sensitivity and specificity parameters were calculated:

	Sensitivity (%)	Specificity (%)
<i>E. coli</i> O157	95.7	97.2
<i>E. coli</i> O145	97.9	100

All samples positive for O157 were negative for O145 and vice versa, demonstrating the serospecificity of the test.

Evaluation of diagnostic performance in early stages of infection

To evaluate the performance of CHEMSTRIP® *E. coli* O157/O145 kit in early stages of the infection, serum samples obtained from patients with HUS or ABD and positive isolation for O157 or O145 were analyzed, and the percentage of detection was calculated at different intervals of time according to the days elapsed between sampling and the onset of symptoms (days post-onset of symptoms, dpos). When the samples obtained before 3 dpos were analyzed, the detection percentage was 80% and after 3 dpos this percentage rose to 100%, demonstrating the usefulness of the kit for the early diagnosis of *E. coli* O157 or O145 infection.

Matrix effect evaluation

CHEMTEST™

The different types of matrices (whole blood and plasma) and anticoagulants (heparin and sodium citrate) evaluated do not affect the test result relative to serum samples.

For technical assistance

For any questions, please contact tel (+54 11) 5353 6066 or through info@chemtest.net / chemtest.net

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