VETLIS® Aujeszky gE bELISA Porcinos

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VETLIS®
Aujeszky gE bELISA Porcine

Presentations / Components

Components	Description	15-VL07-2P 2 microplates	15-VL07-5P 5 microplates
Microplate	96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigen, sealed and stored in dry	2 microplates (192 tests)	5 microplates (480 tests)
HRP anti-gE conjugate (blue)	Ready to use (horseradish peroxidase conjugated anti-gE antibodies)	1 x 30 ml	1 x 62.5 ml
Substrate-chromogen solution	Ready to use [3,3',5,5'-Tetramethylbencidine (H ₂ O ₂)] – STORE IN THE DARK	1 x 30 ml	1 x 60 ml
Stop solution	Ready to use (contains HCl 2%) - CORROSIVE	1 x 20 ml	1 x 40 ml
Wash solution	10X concentrate	1 x 60 ml, 10X (c. s. p. 0.6 litros)	1 x 120, 10X (c. s. p. 1.2 litros)
Sample diluent (violet)	Ready to use - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 40 ml	1 x 100 ml
Positive control (red)	CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 1 ml	1 x 2 ml
Negative control (green)	CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 1 ml	1 x 2 ml
Instruction manual		1	1

Name and application

VETLIS®

Aujeszky gE bELISA Porcine

Blocking ELISA kit for the detection of anti-gE antibodies to Aujeszky's disease virus in porcine serum.

The VETLIS® Aujeszky gE bELISA Porcine kit is a blocking ELISA kit (Enzyme Linked Immunosorbent Assay) for the detection of specific antibodies against the gE protein of the Aujeszky's disease virus in porcine serum.

Introduction

Aujeszky's disease (pseudorabies) is a highly contagious, economically significant disease of pigs. This viral infection tends to cause central nervous system (CNS) signs in young animals, respiratory illness in older pigs, and reproductive losses in sows. Aujeszky's disease results from infection by Suid Herpesvirus 1 (SuHV-1), which is also known informally as Aujeszky's disease virus (ADV). The VETLIS® Aujeszky gE bELISA Porcine kit is designed to detect the presence of anti-ADV gE/gI antibodies in swine serum. It provides a rapid, simple and specific method for differentiating ADV infected animals from those vaccinated with virus strains deleted in the gene that express gE/gI. Animals infected with ADV produce antibodies against ADV gE/gI proteins. However, swine vaccinated with gE/gI deleted vaccines do not produce antibodies against gE/gI proteins. The presence of gE/gI antibodies are evidence of viral infection.

Principle of the technique

VETLIS® Aujeszky gE bELISA Porcine is based on the use of ELISA plates coated with ADV antigens and an indicator system to detect the presence or absence of anti-gE antibodies in serum samples using a monoclonal anti-gE

antibody conjugated to horseradish peroxidase (HRP).

The first step involves performing a 1:2 dilution of the serum on the plate coated with the ADV antigens. During the first incubation, anti-ADV antibodies present in the serum, including those produced against gE, react with the antigens on the plate. After washing, an HRP (horseradish peroxidase)-labeled monoclonal anti-gE antibody conjugate is added to the well and allowed to compete for the gE viral antigen during a second incubation. If anti-gE antibodies are absent in the serum, the gE antigen is available to react with the conjugate. Conversely, if anti-gE antibodies are present in the serum, conjugate binding is blocked. After this incubation period, the unreacted conjugate is removed by washing, and a substrate/chromogen solution is added. In the presence of the enzyme, the substrate is converted into a product that reacts with the chromogen, generating a blue color. The reaction is stopped by adding the stopping solution, and a color change to yellow is observed. The results are read using a microplate spectrophotometer by measuring the absorbance at 450 nm. The measured Abs is inversely proportional to the amount of specific anti-gE antibodies present in the serum sample.

Storage and expiration

Store: between 2 and 8°C.

Transport temperature: 4 to 15°C for 72 hours. Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 18 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

- High precision multichannel pipette and/or micropipettes (50, 100 and 200 µl).
- Containers for multichannel pipette.
- Disposable pipette tips.

- Microplate reader for absorbance determination at 450 nm.
- Timer.
- Tubes for the dilutions of the controls and samples.
- Container/burette for the preparation of the wash solution.
- Distilled (Type II) or ultrapure (Type I) water or any similar high quality water.
- Microplate covers (lid, aluminium foil or adhesive).

Precautions

1. Carefully read and strictly follow all instructions detailed in the instruction manual included in the kit.
2. Store the kit and all reagents at 2-8 °C.
3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use. The required amount of Substrate-chromogen solution and Conjugate should be removed from the refrigerator only 30 minutes before use (10 ml/plate).
4. Handle all materials and reagents according to the Good Laboratory Practices.
5. Serum samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.
6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.
7. Do not mix components or instruction manuals from different test kits batches.
8. Care should be taken to avoid contamination of kit components.
9. If multichannel pipette containers are reused, they should be washed thor-

oughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls in duplicate in each test.

Reagents preparation

a. Wash solution 1X

Prepare 150 ml of 1X Wash Solution (enough for a microplate) by adding 15 ml of 10X Wash Solution to 135 ml of distilled (Type II) or ultrapure (Type I) water. Mix well before using. Once prepared, the 1X Washing Solution can be stored at 2-8°C for 7 days.

Note: if crystals appear in the 10X Wash solution, the bottle must be heated not above 45°C and shake well until completely dissolve.

It is recommended to prepare sufficient quantity to process the plates or wells of the day.

b. Preparation of controls

-Positive control (PC): in tubes of 1.5 ml add 100 µl of the PC to 100 µl of Sample diluent (1/2 dilution). Mix well before using. In each microplate the PC should be added at least in duplicate (100 µl per well).

-Negative control (NC): in tubes of 1.5 ml add 100 µl of the NC to 100 µl of Sample diluent (1/2 dilution). Mix well before using. In each microplate the NC should be added at least in duplicate (100 µl per well).

c. Sample preparation

Serum samples should be diluted 1/2 using the Sample diluent. Add 100 µl of serum to 100 µl of Sample diluent. Use 100 µl of diluted sample per well.

Alternatively, perform the 1/2 dilution of the controls and samples directly on the plate: add 50 µL of Sample diluent and then 50 µL of the control or sample. Homogenize the mixture by pipetting up

and down at least four times.

Notes:

-Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bring to room temperature and homogenize before performing the analysis.

-Samples can be analyzed in a single well. However, it is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

Test procedure

1. Remove the plate from its pouch.

Important: do not remove the plate from its pouch until it has reached room temperature. If the entire plate is not used, separate only the 8-well strips needed to analyze the controls and samples. Store the rest of the strips, along with the desiccants, in the original foil bag hermetically sealed and store in the refrigerator (2-8°C).

2. Add 100 µl per well of each serum sample. See Item Sample preparation. Be careful not to spill samples from well to well.

Note: samples can be processed in a single well or in duplicate. However, to obtain more reliable results it is recommended to run samples at least in duplicates.

3. Add in duplicate 100 µl per well of the NC and PC. See Item Preparation of controls.

4. Cover the plate and incubate at $37 \pm 1^\circ\text{C}$ for 60 minutes.

Note: cover the plate with a plastic lid or aluminum foil to prevent evaporation.

5. Wash/rinse the plates 3 times with 200 µl of 1X Wash Solution per well. See item Reagents preparation.

Note: remove the liquid completely after each washing by dump and tapping onto absorbent material. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well.

7. Cover the plate and incubate at $37 \pm 1^\circ\text{C}$ for 30 minutes.

Note: cover the plate with a plastic lid or aluminum foil to prevent evaporation.

num foil to prevent evaporation.

8. Repeat the washes as indicated in step 5.

9. Add 100 µl of Substrate-chromogen solution to each well, and incubated for 10 minutes (± 1 min) at room temperature ($25 \pm 1^\circ\text{C}$) and in the dark (for example, cover the plate with aluminum foil). Shake manually. Start counting the time after filling the first well.

Note: substrate-chromogen solution should be warmed up for 30 minutes at room temperature (20-25°C) before use (10 ml/plate).

10. Stop the reaction by adding 50 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. Be careful not to spill. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 9).

11. Measure the absorbance at 450 nm (Abs450) immediately using a plate spectrophotometer. There is no need to perform a blank reading.

Important: respect incubation times.

Calculations

Expression of results

Results are expressed by calculating the S/N ratio (Abs450 of the sample over the Abs450 of the CN). The S/N ratio is calculated using the following formula:

$$S/N = \frac{\text{Sample Abs450} *}{\text{NC mean Abs450}}$$

* If the sample was run in duplicate, consider the average of the Abs450 value.

Validity criteria

To confirm the validity of the test the following criteria must be met:

- The Abs450 values of the duplicate controls (PC and NC) and the samples must not differ by more than 20% from the average of the corresponding duplicate.
- The mean of the S/N ratio of the PC must

be less than 0.3 (S/N of PC < 0.3).

- The mean Abs450 of the NC must be greater than 0.5 (Abs450 of NC > 0.5).

If any of these criteria are not met, the test must be considered invalid. In this case, the test must be repeated, strictly following the procedure detailed in the kit's instruction manual.

Interpretation of results

The test sample results should be interpreted as follows:

Serum (dilution 1/2)

S/N	Interpretation
$S/N \leq 0.65$	Positive
$S/N > 0.65$	Negative

Example of result calculation and interpretation

1) Abs450 of PC (duplicate): 0.301 and 0.305

Mean Abs450 = $(0.301 + 0.305) / 2 = 0.303$

2) Abs450 of NC (duplicate): 1.521 and 1.501

Mean Abs450 = $(1.521 + 1.501) / 2 = 1.511$ (valid)

Mean S/N = $0.303 / 1.511 = 0.20$ (valid)

3) Abs450 of Sample: 0.431

S/N value of the sample = $(0.431 / 1.511) = 0.285$

4) Result interpretation: Positive

To facilitate the calculations, the evaluation of the validity criteria of the test and the interpretation of the results, a calculation sheet is available on the company's website (<https://www.chemtest.net/>)

elisa-aujeszky.php) or request it via email to info@chemtest.net.

Important: always check that the version number of the instruction manual provided in the kit matches the version of the spreadsheet.

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Chemtest Argentina S.A.